


Order code: DB686 Nano

# DB686 Nano LED curing light

Small, just right;  
Use, just handy.

Small body, powerful function  
Light curing & caries detection

# 2 in 1



### Technical:

Adapter: AC100-240V, DC5V/1.5A, 50/60Hz

Light curing mode:

Wavelength: 380-520nm

Curing time: 10s/20s

Curable depth: 8mm


Caries detection mode:

Wavelength: 380-420nm


LED light  
measure system




Ultra-small all-metal body  
outer diameter of only 13.5mm




8 mm deep curing depth  
Curing effect into every corner



10mm diameter  
large area light emission




Wavelength: 380-520nm  
Suitable for most of the resins  
in the market




### Detection mode:

Simple, Effctive  
Innovative caries detection



Green:Healthy



Red:Bacterial dentin

Nano can maintain strength even if  
it is not focused on the center.  
Curing at any angle can also ensure the effect

